

# POWERED LAPAROSCOPIC LINEAR STAPLER

# Smart, Precise And Reliable Stapling With One Hand



## **PERFORMANCE**

#### **Consistent Staple line**

 Power-controlled stapling delivers consistent staple formation across broad range of tissue thickness

#### **Continuous Powered Articulation**

Up to 45 degree in each direction to address complicated surgery

# **MANEUVERABILITY**

#### **Precise Control**

- Powered articulation, clamp, firing and release
- Tissue adaptive technology, allows smart control on compression force, time, and firing speed

### **Smart Power Management**

- Smart battery management for long procedures
- Post-op automatic battery drain for better safety
- Provide both non-rechargeable and rechargeable batteries for environment protection

# **SAFETY**

#### **Dry Firing Prevention**

 Stapler firing is not allowed in case of no reload or used reload in the cartridge to avoid dry firing and possible accident

#### **Tissue Adaptive Technology**

 Stapler firing is not allowed when the tissue thickness exceeds limitation to avoid abnormal staple formation

#### **Emergency Gears for Articulation and Firing**

Easy handling in case of emergency to ensure patient safety

# **SIMPLICITY**

#### **Easy Operation**


Ergonomic design of handle and button for easy operation

#### Single Hand Operation

Easy operation with one hand



# **Linear Cutter Model**



# **Linear Cutter Reloads Model**



Reload Code	Staple Line Length	Staple Height	Closed Staple Height	Cartridge Color	Number of Staples	
OCR4520	45 mm	2.0 mm	0.75 mm	Grey	66	
OCR4525	45 mm	2.5 mm	1.0 mm	White	66	
OCR4530	45 mm	3.0 mm	1.3 mm	Bronze	66	
OCR4535	45 mm	3.5 mm	1.5 mm	Blue	66	
OCR4538	45 mm	3.8 mm	1.8 mm	Gold	66	
OCR4540	45 mm	4.0 mm	2.0 mm	Green	66	
OCR6025	60 mm	2.5 mm	1.0 mm	White	90	
OCR6030	60 mm	3.0 mm	1.3 mm	Bronze	90	
OCR6035	60 mm	3.5 mm	1.5 mm	Blue	90	
OCR6038	60 mm	3.8 mm	1.8 mm	Gold	90	
OCR6040	60 mm	4.0 mm	2.0 mm	Green	90	

Surgnova Healthcare Technologies

Tel: +86 10-6411-8265 Fax: +86 10-6411-3261 E-mail: mkt@surgnova.com Web: www.surgnova.com





